



How to Use the EpiPen

How To Use EpiPen®

The EpiPen® Auto-Injector is a disposable, pre-filled automatic injection device that administers epinephrine in the event of a severe allergic reaction.

Remove the EpiPen® Auto-Injector from the carrier tube and follow these 2 simple steps:



- Hold firmly with orange tip pointing downward.
- Remove blue safety cap by pulling straight up. Do not bend or twist.



- Swing and push orange tip firmly into mid-outer thigh until you hear a 'click'.
- Hold on thigh for several seconds.



Built-in needle protection.

After injection, the orange cover automatically extends to ensure the needle is never exposed.

Call 911 after using EpiPen®.

<https://www.epipen.ca/en/about-epipen/how-to-use>